

Blæreskanner med probe som tåler fall fra 1,2 meter. Velges av alle som vet hva det koster å miste en probe i gulvet. Spør oss om pris.

Produktnummer: CM80600

PBS V5.1 er en ny blærescanner med garantert høy nøyaktighet, samt muligheten til å velge mellom enkel og avansert ekspert måling av urinblæren. Den enkle målingen gir intuitive fargekoder som viser når blæren er funnet og målingen kan starte. Ingen behov for opplæring. Ekspert modus gir et langt bedre bilde enn hva som er vanlig og til og med urinstener kan oppdages. OBS!! Proben ødelegges IKKE ved at den mistes i gulvet (1,2 meter). Innebygd printer. Norsk tekst. CE, RoHS, samsvarserklæring og Norsk hurtig manual. Stor skjerm som kan vippes opp på 8,4". Innebygd printer. Uten vifte og er lydløs.



Transducer frequency	2.6MHz
Bladder Volume Range	0-999ml
Scan Mode	Male, Female, Pediatric
Realtime scanning and aiming	Yes
Automatically Outline bladder edge	Yes
Scan Planes	12/24 planes
Scan depth	>18 cm
Sway Angle	95°
Accuracy	±10%,±10ml(0~150ml); ±10%(150~999ml)
Data transmission	USB
Power supply	Lithium Battery
Standby time	4h
Display screen	8.4 TFT
Printer	internal
Device Dimension	300 x 250 x 80mm
Net weigh	1.8kg
Scan time	5s
Working temperature	5-40

Les mer om produktet her:

<https://www.caremed.no/product?number=CM80600>